

# SpiralTOF-TOF

## Synthetic Polymer Structure Analysis

### Poly Methyl Methacrylate (PMMA)

#### Introduction:

The JMS-S3000 “SpiralTOF™” is a MALDI-TOFMS that incorporates an innovative SpiralTOF ion optics system. This system is available with a TOF-TOF option that can acquire high-energy collision-induced dissociation (CID) product ion spectra for monoisotopically selected precursor ions.

In this work, we analyzed Poly Methyl Methacrylate (PMMA) shown in Fig. 1 by using the JMS-S3000 SpiralTOF with the TOF-TOF option. The resulting high-energy CID data was then processed using the Polymerix™ (Sierra Analytics, Inc., [http://masspec.com/](http://massspec.com/)) analysis software.

#### Samples:

Polymer: PMMA  
 Matrix agent: 2,5-Dihydroxybenzoic acid (DHB)  
 Cationization agent: NaI

#### Results and Discussion:

The MALDI mass spectrum of PMMA and the product ion spectra for  $m/z$  1525.8 ( $n=15$ ,  $[M+Na]^+$ ) and  $m/z$  2326.2 ( $n=23$ ,  $[M+Na]^+$ ) are shown in Fig. 2. These product ion spectra show the full range of ions from  $m/z$  23.0 for  $[Na]^+$  to the precursor ion  $m/z$  1027.7 and  $m/z$  2326.2, respectively. Also worth noting, these product ions are all monoisotopic because the precursor ion was monoisotopically selected, which greatly simplified the resulting MS-MS spectrum.




Figure 1. Structural formula of PMMA.




Figure 2. MALDI mass spectrum of a) PMMA, b) product ion spectrum of  $m/z$  1525.8 ( $n=15$ ,  $[M+Na]^+$ ), and c) product ion spectrum of  $m/z$  2326.2 ( $n=23$ ,  $[M+Na]^+$ ).

The enlarged region ( $m/z$  730-850) in Fig. 2b shows that there are at least four different product ions present that could result from the monomer repeat unit (100u,  $C_5H_8O_2$ ). Therefore, this high-energy CID data for PMMA suggests that there are four possible fragmentation pathways.

Based on previously published work<sup>1</sup>, we hypothesized that the structural formulas for these different product ion series were most likely the structures shown in Fig. 3a-d. These formulas were then used in the Polymerix™ software which resulted in the product ion series assignments shown in Fig. 5. These results supported our hypothesis for each structure and showed that the JMS-S3000 high-energy CID data fully represents the structural fragmentation expected for PMMA.

### Conclusions:

Structural analysis of synthetic polymers such as PMMA can easily be done by using the JMS-S3000 TOF-TOF method (high-energy CID, monoisotopic precursor selection) with the Polymerix™ analysis software.

### Reference:

(1) Liang Li, "MALDI Mass Spectrometry for Synthetic Polymer Analysis", John Wiley & Sons, Inc., United States of America (2010).




Figure 3. Structural formula of product ions series.




Figure 4. Polymerix analysis results of  $m/z$  1525.8 ( $n=15$ ,  $[M+Na]^+$ ) (upper) and  $m/z$  2326.2 ( $n=23$ ,  $[M+Na]^+$ ) (lower).