

Data Integrity Compliance of JEOL NMR

JEOL NMR systems have had all FDA 21 CFR Part 11 support functions available in the standard Delta software since version 4.3 was released in 2004. All newer versions of the software have maintained and added to this functionality.


An essence of GMP management is regularization and documentation, and thus the software is required to provide 4 basic functions as detailed below.



In a JNM-ECZ series NMR system, a built-in spectrometer control computer in the spectrometer and an operation PC are completely independent of each other. This system architecture allows for very secure operation of data management. An NMR data (JDF format) file generated by an JNM-ECZ series spectrometer includes not only the traditional NMR data but also all kinds of metadata, e.g. measurement parameters, data processing parameters, audit trails, and electronic signature(s) in single file. It is not necessary to connect multiple files, making data management more straightforward.

ALCOA-plus principles support

Modern data integrity guidance requires 5 principles (ALCOA) as well as 4 additional principles (CCEA, also known as ALCOA-plus) to be followed by the system.



ALCOA	Attributable	Delta software automatically records when/who/where/what have been done during all data collection and all data processing steps.
	Legible	Delta software equips File Information Viewer to browse records and metadata.
	Contemporaneous	Delta software records data collection and data processing with time stamps at the same time of operation.
	Original	Delta software does not allow alteration of data without records nor any correction of records.
	Accurate	Delta software automatically records time stamp, user name, computer name, IP address and operational details, and thus there is no way to tamper with the data or falsify information.
CCEA	Completeness	NMR data file of Delta software, which is binary data, includes all metadata including pulse program, measurement parameters, data processing details, and etc.
	Consistent	Delta software records all data collection including interrupted experiment, and also records all data processing details including reprocessing.
	Enduring	Remote Maintenance Tool, which is provided with Delta software, equips a function to back up all data regularly.
	Available	Delta software is forwards compatible in principle, and compatibility of the file information browsing feature is always kept.

Certain products in this brochure are controlled under the "Foreign Exchange and Foreign Trade Law" of Japan in compliance with international security export control. JEOL Ltd. must provide the Japanese Government with "End-user's Statement of Assurance" and "End-use Certificate" in order to obtain the export license needed for export from Japan. If the product to be exported is in this category, the end user will be asked to fill in these certificate forms.

